

UBE2K-ubiquitin-conjugating enzyme E2K (UBC1 homolog, yeast) Gene

Reference: PTXBC022804

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E2-25K; HIP2; HYPG; LIG; UBC1; ubiquitin-conjugating enzyme E2 K; E2 ubiquitin-conjugating enzyme K; E2(25K); HIP-2; huntingtin-interacting protein 2; ubiquitin carrier protein; ubiquitin conjugating enzyme E2K; ubiquitin-conjugating enzyme E2(25K); ubiqu		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2K (UBC1 homolog, yeast)
GENE ID	3093
ORIGIN SPECIES	Human
ACCESSIONS	BC022804
VECTOR	pDONR223
NCBI SYMBOL	UBE2K
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaacatcgcggtgcagcgaatcaagcgggaggtcaaggagtgctgaagagcgaggagacgagcaaaaaatcaataaagtagatctttagatgagaattttacagaattaaggaggagaatagcaggacctccagacacaccatatgaaggagggaagataccaactagaga
taaaaataccagaaacataccatftaatccccctaaaggtccgggtttatcactaaaatatggcatcctaataattagttccgtcacaggggctatttggatctctgaaagatcaatgggcagctgcaatgactctccgcacggtatttgcattgcaagcactattggcagctgcagagcca
gatgatccacaggatgctgtagtagcaaatcagtacaacaaaatccgaaatgttcaaacagacagctcgacttgggcacatgtgtatgctggagcaccagtttctagtcagaatacaccaaaaaatagaaaacctatgtgctatgggctttagatgaatgcagtaatagtgcccttgc
ttcaaatcatgggatgtagagactgcacagaattgctctgagtaactga