

TMEM111-transmembrane protein 111 Gene

Reference: PTXBC022807

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM111; POB; ER membrane protein complex subunit 3; 30 kDa protein; partial optokinetic response b; transmembrane protein 111		

For research use only.

GENE DESCRIPTION	transmembrane protein 111
GENE ID	55831
ORIGIN SPECIES	Human
ACCESSIONS	BC022807
VECTOR	pDONR223
NCBI SYMBOL	TMEM111
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggccagaaactgttctcgaactccaacatccgctctgggtgtgttaccatcgttatcatcactttctcgttagcatgatccgactacgtgtccatcctgctgcagagcgacaagaagtcacccaggacaagaatctgcagtcgaagtccttaattcgaagcagagtcctca
gggaaatggaaatatactcccaacagcttcttgacacgaaatattttcaacaacccagagatggattttcaaaaaaactaaacggaagtagtgcacacttctcctatgactgactcctactatgttgacagacatgatgaaagggaatgaacaaatgtcctccctatgattctattg
gtggatgatcaacatgacattctcaggctttgtcacaaccaaggtccatttccactgacctcgggttttaagcctatgttacagcaaggaatcgagctactacattagatgatcctgggtgagttctgcatcctggtacttctcctcaatgtatttgggcttcggagcatttactctgattctggg
ccaagataatgccgctgaccaatcacgaatgatgcaggagcagatgacgggagcagccatggccatgccgcagacacaaacaaagcttcaagacagagtggaagcttggagctgacggatcaccagtgggcactagatgatgtcgaagaagagctcatggccaaagacctcca
cttcgaaggcatgttcaaaaaggaaatfacagacctatattttga