

ZNF32-zinc finger protein 32 Gene

Reference: PTXBC022842

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KOX30; zinc finger protein 32; C2H2-546; zinc finger protein 32 (KOX 30); zinc finger protein KOX30		

For research use only.

GENE DESCRIPTION	zinc finger protein 32
GENE ID	7580
ORIGIN SPECIES	Human
ACCESSIONS	BC022842
VECTOR	pDONR223
NCBI SYMBOL	ZNF32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttgattccaacagctaccctgtggactgtcatggaagatagccagaatgtagcgtttcattatgtgactgaagccaccacaaatataccactctgaggctacagatcctcaagctgggatatccaaattcttcagaagagagaagctggaacaaaatccccagattc
gaagacactacaggaagattcaccctggagtgagacaaaggtctatgagtgccaggagtggtgaaatccttcggcaaaaaggtagttaacgttacatgagagaatccacactggtcaaaagccttttgagtgacccactgtgaaaaagcttcaggccaaaggcaatctgttaca
catcaacggatacacacgggagagaagccttaacgtgcaaggtgtggtgaaagcttcagtaacgaggtagtcgtgtccacgagagactccacactggacagaaacctacgagtggtctattgtcagagaagcttcaggaaatcagagtaacctgtgttcacaggagagttc
acagtggtgagaagcctatagatgtgatcagtggtgaaagccttcagtcagaaaggaagccttaattgttcacatcagagtcacacagggcctgaagccctatgcctgtaccagtgagggaagagttccaccagggggaattgtatttcgatggcaaaatccacacagagagac
accctatctgtgcccagtggtgaaagccttcaccagagaggagtcgtgtgcaccagcgaagctgtcacagaggctcaccccttga