

THAP6-THAP domain containing 6 Gene

Reference: PTXBC022989

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	THAP domain-containing protein 6; THAP domain containing 6		
For research use only.			

GENE DESCRIPTION	THAP domain containing 6
GENE ID	152815
ORIGIN SPECIES	Human
ACCESSIONS	BC022989
VECTOR	pDONR223
NCBI SYMBOL	THAP6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaaatgctcgcgcaatggatgtgcttgccttgccaaatcgaagttaaaaggactgacatttcacgtattccccacagatgaaacatcaaaaggaaatgggtattagcaatgaaaagacttgatgtgaatgcagccggcattggggagcctaaaaaggagatgtgtgtg
ttcgaggcactttaagaagacagattttgacagaagtgctccaaatattaaactgaaacctggagtcataccttctatctttgattctccatcacctacaggggaaaagagaaaaacttcattgtagaaaaaacttcacctcaaaaccgttcagccactaactacaatcaccatctgtgtgtg
cttcctcatgtattgaagaattccaatccagttcatgtttgaacatagctacagtgtgaatggacagtcctcaaaagaaacttaagcataaattagatcatgtgatcgccgagctagaggatatacaaaaggaaagtctacggaatgttttagaccgagaaaaacgttttcagaaatcattgaggaagacaat
cagggaaattaaaggatgaatgtctgacagccaagaacagcaaatagacttgacactttctgttgggactgtgtcaggagagcatagaacaggactatattcatga