



Reference: PTXBC023152

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GN-2; GN2; glycogenin-2; glycogenin glucosyltransferase; glycogenin 2		
For research use only.			

GENE DESCRIPTION	glycogenin 2
GENE ID	8908
ORIGIN SPECIES	Human
ACCESSIONS	BC023152
VECTOR	pDONR223
NCBI SYMBOL	GYG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]