



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

VAMP1-vesicle-associated membrane protein 1 (synaptobrevin 1) Gene

Reference: PTXBC023286

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPAX1; SYB1; VAMP-1; vesicle-associated membrane protein 1; vesicle-associated membrane protein 1 (synaptobrevin 1); vesicle associated membrane protein 1		

For research use only.

GENE DESCRIPTION	vesicle-associated membrane protein 1 (synaptobrevin 1)
GENE ID	6843
ORIGIN SPECIES	Human
ACCESSIONS	BC023286
VECTOR	pDONR223
NCBI SYMBOL	VAMP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgctgcctgacgcacctgctgctgaagggacagagaaggggactgcccgactggggggctcccttgcctctctactcaatgacgacgactaacagcagactacagcaaacccaggcacaaggctggaggaggtggtggacatcatcgtgtgaactgggacaaggtcctggagaggggacagaagctgtcagactgctgatgaccagagctgactgcttgcagcagaggacatcacaaatttgagagcagctgctgcaaaagttaagagaggaggtattgtgtggaaaaactgcagaatgatgatcatcgtctggagcgaatcgtgccaatcatcgtgtgtagttattgaagtaagatcctgta