



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ST3GAL6-ST3 beta-galactoside alpha-2,3-sialyltransferase 6 Gene

Reference: PTXBC023312

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SIAT10; ST3GALVI; type 2 lactosamine alpha-2,3-sialyltransferase; CMP-NeuAc:beta-galactoside alpha-2,3-sialyltransferase VI; alpha2,3-sialyltransferase ST3Gal VI; sialyltransferase 10 (alpha-2,3-sialyltransferase VI); ST3 beta-galactoside alpha-2,3-sialyl		

For research use only.

GENE DESCRIPTION	ST3 beta-galactoside alpha-2,3-sialyltransferase 6
GENE ID	10402
ORIGIN SPECIES	Human
ACCESSIONS	BC023312
VECTOR	pDONR223
NCBI SYMBOL	ST3GAL6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr