

GNLY-granulysin Gene

Reference: PTXBC023576

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D2S69E; LAG-2; LAG2; NKG5; TLA519; T-cell activation protein 519; T-lymphocyte activation gene 519; lymphocyte-activation gene 2; lymphokine LAG-2		

For research use only.

GENE DESCRIPTION	granulysin
GENE ID	10578
ORIGIN SPECIES	Human
ACCESSIONS	BC023576
VECTOR	pDONR223
NCBI SYMBOL	GNLY
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctacctgggcccctctgctcctgagccatgctcctgggcaaccaggtctgttctctctctgtgagccctgagctactacgacctggcaagagccacctgctgatgaggagaaatcctgccctgacctggccaggaggccccaggggtgacctgtgacccaaacacag
gagctgggcccgtgactacaggacctgtctgacgatagtccaaaactgaagaagatgggtgataagccaccagagaagtgttccaatgctgcacccgggtgtgttaggacggggaggtcagatggcgcgacgtctgcagaaattcatgaggaggtatcagctagagttaccca
gggctctgtggccgagaaactgccagcagatctgtgaggacctcaggtgtgtatcctctacaggtccctctga