



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TREX1-three prime repair exonuclease 1 Gene

Reference: PTXBC023630

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	3'-5' exonuclease TREX1; AGS1; CRV; DRN3; HERNS; three-prime repair exonuclease 1; 3' repair exonuclease 1; DNase III; deoxyribonuclease III; three prime repair exonuclease 1		

For research use only.

GENE DESCRIPTION	three prime repair exonuclease 1
GENE ID	11277
ORIGIN SPECIES	Human
ACCESSIONS	BC023630
VECTOR	pDONR223
NCBI SYMBOL	TREX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr