



Reference: PTXBC023646

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein spinster homolog 3; SPNS sphingolipid transporter 3 (putative); spinster homolog 3; sphingolipid transporter 3 (putative)		

For research use only.

GENE DESCRIPTION	spinster homolog 3 (Drosophila)
GENE ID	201305
ORIGIN SPECIES	Human
ACCESSIONS	BC023646
VECTOR	pDONR223
NCBI SYMBOL	SPNS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]