

IGLV6-57-immunoglobulin lambda variable 6-57 Gene

Reference: PTXBC023973

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	immunoglobulin lambda variable 6-57
GENE ID	28778
ORIGIN SPECIES	Human
ACCESSIONS	BC023973
VECTOR	pDONR223
NCBI SYMBOL	IGLV6-57
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctgggtccactacttctcaccctctcgtcactgcacaggttcttgggccaattttatgctgactcagccgcaactgtgtcggagtcctcggggaagacggttaacctctctgcacccgcagcagtgccagcattgccagcaactatgtcagtggtaccagcagcgggg
ccgtgccccaccactgtgatctatgaggataaccaagacccctcgggtccctgacgggttctctggtccatcgacagctcctccaactctgcctccctcaccatctctggactgaagactgaggacgaggctgactactactgctcgtcttatgatagcagcaatcacacagtgtctcag
acccatggggaagtgagacagaaactccccagagcatctctactcgggccagttcagcctgtctccaccggagagggttag