



Reference: PTXBC024043

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAN2 poly(A) specific ribonuclease subunit; PAN2 polyA specific ribonuclease subunit; PAN2 poly(A) specific ribonuclease subunit homolog; PAB-dependent poly(A)-specific ribonuclease subunit PAN2; USP52; PAB-dependent poly(A)-specific ribonuclease subunit		

For research use only.

GENE DESCRIPTION	PAN2 polyA specific ribonuclease subunit homolog (S. cerevisiae)
GENE ID	9924
ORIGIN SPECIES	Human
ACCESSIONS	BC024043
VECTOR	pDONR223
NCBI SYMBOL	PAN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]