

C22orf32-chromosome 22 open reading frame 32 Gene

Reference: PTXBC024237

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0466 protein C22orf32, mitochondrial; C22orf32; DDDD; EMRE; dJ186O1.1; essential MCU regulator, mitochondrial; single-pass membrane protein with aspartate-rich tail 1, mitochondrial; single-pass membrane protein with aspartate rich tail 1		

For research use only.

GENE DESCRIPTION	chromosome 22 open reading frame 32
GENE ID	91689
ORIGIN SPECIES	Human
ACCESSIONS	BC024237
VECTOR	pDONR223
NCBI SYMBOL	C22orf32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtccgagcggcgtcgtgctagtattggcaccgctcaggtccgggctctccggagcgggcctagcttgaggaaagatggcgatgtctccgccgatggagcggcctcagggcggagcctggtagcgtcggggtcagtcagttaccgcagcggcgccattttgcccaaac
cgggtgaaaatgtcttccttcgctcgtgtgttctccattgtgatccctttctctatgtcgggacactcattagcaagaacttgcctcctacttgaggacatgacattttgtccagaggatgatgatgatgactaa