

CDO1-cysteine dioxygenase, type I Gene

Reference: PTXBC024241

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDO-I; cysteine dioxygenase type 1; cysteine dioxygenase, type I		
For research use only.			

GENE DESCRIPTION	cysteine dioxygenase, type I
GENE ID	1036
ORIGIN SPECIES	Human
ACCESSIONS	BC024241
VECTOR	pDONR223
NCBI SYMBOL	CDO1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacagacgaagtgcgaagccacggaccctggctgatctgatccgcatcctgcaccagctcttggcggcgatgaggtcaatgtagaggaggtgcaggccatcatggaagcctacgagagcgacccaccgagtgaggcaatgtacgccaagttcgaccagtacaggtataccg
gaaatctgtggatcaaggaaatggaaaattaatctgatattctctgttgggtgaaggacatggcagcagttatcatgatcatacctaactcccactgcttctgaagatctacagggaaatctaaaggagacattattgcctggcctgacaaaaatccaatgagatgggtcaagaagtctg
aaagagctctgagggaaaaccagtgtgcctacatcaatgattccgttggcttaccatcgagtagagaacatcagccatacggaaacctgctgtgagccttactgtacagtcacacttttgatcatgccatgcctttgatcaagaacaggacataaaaacaagtacaaatgacattccatagta
aatttggaaatcagaactccaaatgcaactcgggctcgtggagaacaactaa