

IL1F5-interleukin 1 family, member 5 (delta) Gene

Reference: PTXBC024747

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL1F5 (Canonical product IL-1F5a); IL1F5; FIL1; FIL1(Delta); FIL1D; IL-36Ra; IL1HY1; IL1L1; IL1RP3; IL36RA; PSORP; PSORS14; interleukin-36 receptor antagonist protein; IL-1 related protein 3; IL-1F5 (IL-1HY1, FIL1-delta, IL-1RP3, IL-1L1, IL-1-delta); IL-1		

For research use only.

GENE DESCRIPTION	interleukin 1 family, member 5 (delta)
GENE ID	26525
ORIGIN SPECIES	Human
ACCESSIONS	BC024747
VECTOR	pDONR223
NCBI SYMBOL	IL1F5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcctgagtgaggcgctgtgctccgaatgaaggactggcattgaaggtgctttatctgcataataaccagcttctagctggaggctgcatgcagggaaggtcattaaaggtgaagagatcagcgtggtcccaatcggtgctggatccagcctgtccccgtcatcctgggtgtccagggtggaagccagtgctgtcatgtgggtgggagagccgacttaacactagagccagtgaaacatcagagctctatcttggtgccaaggaatccaagagcttcaccttaccggcggaatggggtcactccagcttcagtcggctgacctaccgggctgttctgtgcacgggtgctgaagccgatcagcctgtcagactcaccagcttcccgagaatgggtggctggaatgccccatcacagacttctactccagcagtgtagtag