

GAGE5-G antigen 5 Gene

Reference: PTXBC024914

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT4.5; G antigen 5; GAGE-5; cancer/testis antigen 4.5; cancer/testis antigen family 4, member 5		

For research use only.

GENE DESCRIPTION	G antigen 5
GENE ID	2577
ORIGIN SPECIES	Human
ACCESSIONS	BC024914
VECTOR	pDONR223
NCBI SYMBOL	GAGE5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttggcgaggagatcgacattattggcctagaccaaggcgctatgtacagcctctgaagtgttggcctatgcggccgagcagttcagtgatgaagtgggaaccagcaacacctgaagaaggggaaccagcaactcaacgtcaggatcctgcagctgctcaggaggag
aggatgaggagcattcgcaggtcaaggccgaagcctgaagctgtagccaggaagggtcaccacagactgggtgtgagtgtagatgtgctctgatgggcaggagatggacccgcaaatccagaggaggtgaaaacgcctgaagaaggtgaaaagcaatcacagtgtta
a