

FBXO2-F-box protein 2 Gene

Reference: PTXBC025233

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FBG1; FBX2; Fbs1; OCP1; F-box only protein 2; F-box gene 1; organ of Corti protein 1; F-box protein 2		

For research use only.

GENE DESCRIPTION	F-box protein 2
GENE ID	26232
ORIGIN SPECIES	Human
ACCESSIONS	BC025233
VECTOR	pDONR223
NCBI SYMBOL	FBXO2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacggagacgggtgaccagagagcgtggccagcccgaggaggaagcccgaggagcagccagaggagcagtgctgaggagggcggccggcggcggcgtacgtggagagctgcccagcggcgtgctg
cgctgctggccgactgcccggccggagctgggtgagggcctgcccgtgggtgctgctggaagagcgtgggtgagcggcggccgctggtgctgctcaagtgccagcaggaggcgtggtgcccaggggcggcgtggaggaggcgcgaccactggcagcagttcta
cttcctgagcaagcggcggcgaacctctgctgaacctgtgtggggaagaggacttgaaggcgtgggtgacgtggagcatggtgggacggcgtggagggtggaggagctgctggagacagtggggtggagttcaccacgatgagagcgtcaagaagtacttgcctcctctt
gagtggtgtcgaagcacaggtcattgacctgacgtgagggctactgggaggagcgtgctggacacgactcagccggccatcgtggtgaggagcgtgactcgggcccagcgcgacgctggtgctctacgagctcacgttaagctactgtccgagcagagaacgtgctgctg
agttcagcagcggcaggtggcagtgcccaagacagtgacggcggggcgtggatggagatctccacaccttcaccgactacggccggggcgtccgttcgtccgttcgagcacggggggcaggactcgtctactggaaggcgtggtcggggcccggtgaccaacagcag
cgtgtgggtagaacctga