

hCG_2018279-hypothetical protein LOC100127888 Gene

Reference: PTXBC025345

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC100127888		
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC100127888
GENE ID	100127888
ORIGIN SPECIES	Human
ACCESSIONS	BC025345
VECTOR	pDONR223
NCBI SYMBOL	hCG_2018279
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcacactggggcagtgagtgagaagacgagatgccagcaggtgacaacgggtgcagaacaggcagaactgatgaccgggagccccggagcctgagcttgttcacaagaagccgccaggaaggtgccgccaccaccgctggcaccaccagggcagccacacc
ccctctccccaggggcctctgccgagtcctgacggtgtgacaaggacattgtcgaagcctga