

ATP6V1D-ATPase, H⁺ transporting, lysosomal 34kDa, V1 subunit D Gene

Reference: PTXBC025373

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP6M; VATD; VMA8; V-type proton ATPase subunit D; ATPase, H ⁺ transporting lysosomal, member M; ATPase, H ⁺ transporting, lysosomal (vacuolar proton pump); ATPase, H ⁺ transporting, lysosomal 34kDa, V1 subunit D; H(+)-transporting two-sector ATPase, subunit		

For research use only.

GENE DESCRIPTION	ATPase, H ⁺ transporting, lysosomal 34kDa, V1 subunit D
GENE ID	51382
ORIGIN SPECIES	Human
ACCESSIONS	BC025373
VECTOR	pDONR223
NCBI SYMBOL	ATP6V1D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgggcaagaccgaattgaaatcttccctcgcgaatggcagacacatcatgaaggctcgttaaggggagcagacaggtcgaaacctctgaagaaaaatctgatgccttaactcttcgatttcgacagatcctaaagagataatagagactaaaatgttgatggcggaagtg
atgagagaagctgccttttactagctgaagccatgttcacagcaggtgacttcagcactacagttatccaaatgtcaataaagcgcaagtgaagattcgagcgaagaaagataatgtagcaggtgttactttgccagtatttgaacattaccatgaaggaaactgacagttatgaactgactgg
tttagccagagtggtgggaacagttgctaaattaaagaggaattatgccaaagcagtggaactactgttggaactagcttctctgcagacttctttgttacttttgatgaagctattaagataaccaacagggcgtgtaaatgccattgaacatgtcatctccccgattgaacgtactcttgcct
atatcatcacagagctggatgagagagagcgagaagagttctataggttaaagaaaatacagaagaaaaagattctaaaggaaaaatctgagaaggacttggagcaaggagagcagctggagaggtgttgagcctgctaactcttgcctgaagagaaggacgaggtcttctat
ttgaataa