

SOHLH2-spermatogenesis and oogenesis specific basic helix-loop-helix 2 Gene

Reference: PTXBC025383

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SOSF2; SPATA28; TEB1; bHLHe81; spermatogenesis- and oogenesis-specific basic helix-loop-helix-containing protein 2; spermatogenesis associated 28; testicular secretory protein Li 50; testicular secretory protein Li 51; spermatogenesis and oogenesis specif		

For research use only.

GENE DESCRIPTION	spermatogenesis and oogenesis specific basic helix-loop-helix 2
GENE ID	54937
ORIGIN SPECIES	Human
ACCESSIONS	BC025383
VECTOR	pDONR223
NCBI SYMBOL	SOHLH2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttctcaattatctgccagagcactgccagatctcggccaggcacaataagacatcttattgttgagatgtcactgtgggtacctggctgatactgtacagaaactattgcaaacatagcagaagtcaccatcaccatcagtgacacgaaggaggcagcagcgttttggatg
attgcatattcaacatgggttctcttgaagggtgccttctcactaagtcggagagctggaagccatcaagtttaattagatttgcacaaagaaaaatacattcactgtttgttttataatccctgaaaaatttaagggttgatttcagggcattggaatggatattgctttaactgaaccactgacaa
tggaaaaatgagtaattgtgtaaaatactggacaacatgtccctcaaacactgttaagactgaaaacgcaactgggcctgaagaacttgattgcccctgcagaggtctacagcgaacacctgggatatttctactgatctattgcctgctctgaatcttaagggaatggcaatggccttg
aattaaatgcttctgtcagagttcagaaaaacaaaagatctctcttcattcaagcaagaaaactaagaggctgtacaggaagcatagcagcttctgttctgtga