

39515-membrane-associated ring finger (C3HC4) 8 Gene

Reference: PTXBC025394

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	membrane-associated ring finger (C3HC4) 8
GENE ID	220972
ORIGIN SPECIES	Human
ACCESSIONS	BC025394
VECTOR	pDONR223
NCBI SYMBOL	39515
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcatgccactgcatcagatctctgccattccatccagatgccatctctgtagatctacagaagtaagacaaagaaaggagagggaagacagaatgagaagactttggacatttcagtcattcaagcaacatttcaaggctgggagctcctcgtcagcatcagctcc
ggctccgggtgtctctctctctcgcacttcatcagccatccagccagacatctgcaggatctgccactgtgaaggagatgatgagagccctctgatcccccctgccactgcacaggaaagcctccacttctgtgaccaggcctgctgcagcagtgatcaagagctccgacacgcg
ctgctcgcagctctgcaagtatgagttcatcatggagacaaagctgaagccactgagaaaatgggagaagtgtcagatgacgtccagcagcgcaggaaagatcatgtgctcagtgacattccacgtcattgccatcacatgtgtgtctgtgctctgtatgtgctcattgaccgtactgtctgag
gagatcaagcaggggcaggcaacaggaatcctagaatggccctttggactaaattgtgtgtgtgccatcggttcacccggaggactctttttatgtatgttcagtgtaaagtgtatgtgcaattgtgggaagagactcaaggcctataatagatgatctatgttcaaaactgtccagaaaca
agcaaaaagaatattttgaaaaatctccactaacagagcccaacttfgaaaataaacatggacatggaatctgtcattccgacacaaactcttctgttcacagagcctgaagacactggagcagaatactcattcagctctga