

## CD300LG-CD300 molecule-like family member g Gene

Reference: PTXBC025395

|                |                                                                                                                                                                                         |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                     | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                    | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CLM-9; CLM9; NEPMUCIN; TREM-4; TREM4; CMRF35-like molecule 9; CD300 antigen-like family member G; triggering receptor expressed on myeloid cells 4; CD300 molecule like family member g |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | CD300 molecule-like family member g                                                               |
| <b>GENE ID</b>          | 146894                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC025395                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CD300LG                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgcggctctgtgctctgtatgggttgcctgctctccaggttatgaagccctggaggccagaggaaatcagcgggttcgaaggggacactgtgtccctgcagtgcacctacagggagagctgaggaccaccggaagtactggtgcaggaaggtggggtcctctctctg  
ctgctctggcaccatctatgcagaagaagaagccaggagacaatgaaggcaggggtgtccatccgtgacagccgaggaagctctgctcattgtgacccctgtggaacctcaccctgcaagacgctggggagtactggtgtgggtcgaaaaacggggcccgatgagctttactg  
atctctctgttctcttccaggaccctgctgtctccctccctctccacctccagcctctgtgctacaacacgctgagcccaaggcaaaagctcagcaaacccagccccaggtgacttctctgggctctaccggcagccaccacagccaagcaggggaagacaggggctg  
aggccctccattgccagggaactccagtaggggacgaaaggactctcagtagacaggaacctctctcaccagcgacctctctctctgagggagctcccgccccccatgcagctgaactccacctcagcagaggacaccagtcagctctcagcagtgccagctctaaagccc  
agggtgtccatcccgatggtccgatactggccccagtcctggtgctgtgagcctctctcagccgagccctgacgctctctgagccacctgctctgtggagaaaggagctcaacagccacggagacacaggaacgagaagttctgctctcagcgttgcggagga  
aaagggaagccctccagggccctgaggggacgtgactctgagctccctccctccacacatctgaggaggagctgggcttctcgaagttgtctcagcgtag