

CD247-CD247 molecule Gene

Reference: PTXBC025703

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD247 molecule; CD247 antigen, zeta subunit; CD3-ZETA; CD3H; CD3Q; CD3Z; IMD25; T3Z; TCRZ; T-cell surface glycoprotein CD3 zeta chain; CD3Z antigen, zeta polypeptide (TiT3 complex); CD3zeta chain; T-cell antigen receptor complex, zeta subunit of CD3; T-ce		

For research use only.

GENE DESCRIPTION	CD247 molecule
GENE ID	919
ORIGIN SPECIES	Human
ACCESSIONS	BC025703
VECTOR	pDONR223
NCBI SYMBOL	CD247
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagtgaaggcgctttaccgcggccatcctgcaggcacagttgccgattacagaggcacagagctttggcctgctggatcccaactctgctacctgctggatggaatcctctcatctatggtgtcattctactgccttggtcctgagagtgaagttcagcaggagcgacagcccc
ccgcgtaccagcagggccagaaccagctctataacgagctcaatctaggacgaagagaggtacgatgttttgacaagagacgtggccgggacacctgagatgggggaaagccgcagagaaggagaaccctcagggaaggcctgtacaatgaactgcagaaagataagatggc
ggaggcctacagtgaattgggatgaaaggcagcgcggagggggaaggggcacgatggcctttaccagggtctcagtagcaccacagacactacgacgcccttcacatgcaggccctgccccctcgctaa