

## LILRB5-leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 5 Gene

Reference: PTXBC025704

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CD85C; LIR-8; LIR8; leukocyte immunoglobulin-like receptor subfamily B member 5; CD85 antigen-like family member C; leukocyte Ig-like receptor B5; leukocyte immunoglobulin-like receptor 8; leukocyte immunoglobulin-like receptor, subfamily B (with TM and I |                     |                |

For research use only.

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 5          |
| <b>GENE ID</b>          | 10990                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC025704                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | LILRB5                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGCGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaccctaccctcagctcgtgattgctcgggctgagtggtggcccccaggacctgcgtgcaggcagggaccctcccaaacccacctctggctgagccagcctctgtgatagctcggggagcccgtagacctctggtgtcaggggccctggagactgaggagtaccgtct  
ggataaggaggagctccatggcccccgaagagacagaaccactggagcctggagcccaaggcccaagttccacatccatccacgggtgtatgacagtgccagggcgataccgctgctactatgagaccctgcaggctgtgcagagcccagtgacccctggagctgtggtgcagag  
gattctatgcagaaccacactcttttagccctgcccagtcctgtggtggcctcaggaggaatgtgacctccagtggtgatacactggagcagctctcagcttctgttctgtgaggagaagacagaagctcccccaggccatcccaggccctgtt  
ccctgtgggtcccgtagcccccagctgcagggtggaggttcagatgcttactattacaggaaaaaacctcaggtgtggtggaacccagtgacctctggagattctggtccaggcgtgtctaggaagccctccctctgatcccgagggtctgtcgtggcccgaggagcctg  
accctgcagtgctcgtctgatgctggctatgacatattctgtgtacaaaggagggggaacatgacctgtccaggcctctggtccagcagccccaggcctgggtctctccaggccaacttaccctggccctgtgagccgtccacggggccagtagacatgctacggtgcacaaa  
cctctccctaggtgtggtcggcccccagcgacccctggacatcctgatgcaggactgatccctgacatacccgccctctcgggtcagccggcccaagggtggcctcaggagagaaacgtgacctgctgtgtcagtcagatggcatcagatagacacttcttttgaacaggagggggc  
agccatccccctgctgtgtctaaagtcaagtaccagtcttatagacacaggctgaattctccatgagtcctgtgacctcagccagggtggaacctaccgatctacagcgcaatcaggtctaccctacctgctgtccagccctagttaacccagggagctcgtgtcaggaccctc  
tgggatccagcctcctacctacaggctccaccccccacacctggccctgaggaccagccctcaccctcaggggttgatccccagagtggtctgggaaggcaacctgggggtgtgactgggtctcagtgccctcgtcctgtgctgttctcctctctcctcctccgacatc  
ggcatcagagcaaacacagacatcggccatcttctacgctcctgcaggggctcggggccagagcccaaggaccagggcctgcagaagaggccagccaggtgtgacatccaggaggaattctcaatgctcgtggaaggacacacagcccaaggacgggtggagatgg  
atgctcgggtgctgcatctgaagccccccaggatgtgacctacgcccagctacacagcttgacctcagcggggaggaactgagcctctccatccaggaaagggaacctccagctgaaccagcatctacgccccctggccatcactag