



Reference: PTXBC025714

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CALN; CALNA; CALNA1; CCN1; CNA1; PPP2B; serine/threonine-protein phosphatase 2B catalytic subunit alpha isoform; CAM-PRP catalytic subunit; calcineurin A alpha; calmodulin-dependent calcineurin A subunit alpha isoform; protein phosphatase 2B, catalytic su		

GENE DESCRIPTION	protein phosphatase 3 (formerly 2B), catalytic subunit, alpha isoform
GENE ID	5530
ORIGIN SPECIES	Human
ACCESSIONS	BC025714
VECTOR	pDONR223
NCBI SYMBOL	PPP3CA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1