

hCG_1731871-hCG1731871 Gene

Reference: PTXBC025725

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CXorf50B; LINC00246B; NCRNA00246B; family with sequence similarity 226, member B (non-protein coding)		

For research use only.

GENE DESCRIPTION	hCG1731871
GENE ID	653687
ORIGIN SPECIES	Human
ACCESSIONS	BC025725
VECTOR	pDONR223
NCBI SYMBOL	hCG_1731871
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgaagtatgtcatcaaaaggtagcagagcttcctcctgagatttcaagaaagcctctgacctccccgagttagcttgggtctatctgaaggaaactgacccagagcagcactctatgtctgatcagaaaaatgatcctgcctgggttggggcctgacaggcgaccagggcac
accaaagacccggctcctgatgattactctggactcgatctcatgcaggccagctgtgtccccgaggaggtgtctgggaggtgttgagggtgttgaggcaccatttctctaa