

FAM124B-family with sequence similarity 124B Gene

Reference: PTXBC025754

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM124B; family with sequence similarity 124B; family with sequence similarity 124 member B		

For research use only.

GENE DESCRIPTION	family with sequence similarity 124B
GENE ID	79843
ORIGIN SPECIES	Human
ACCESSIONS	BC025754
VECTOR	pDONR223
NCBI SYMBOL	FAM124B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgagacacaggggctctggccatgactgtccatctcttccaactctgggacggctccctctgcagaggactctggaccagctcctggattgcaattgcccagaggtccggctcttcagggtgtctgaacgggccaagctctgtgaaatactgtgaaaagtcacattccaagcgggt
cccgggttccagggatgtccgtgtgtctctctctgcacgaaagccgggagagataggtatcttgcgtctgactctccagcattcgccatggcagtgctacccccacccaggacactcggggaaggctgtgtccctactttttgccaatcaggagttctacagcctggacagtcagct
gcccactctgggggtgagggcaggtgcactgtggctccgagatcctgaggggtgacgtgtactgcagttttgataactatgaagacgcatcagactctacgagatgatcctgcagagagaagcgaccttgcaaaagagcaattttgttctctgtctctatgcctccaagagctttgtctctgc
agctctccctgaagcagctgccccggggaatgtcagtggaacccaaagagttctgggtgctgcagtttaagggttcaagagatcgccagtttagtgcctctgctaccaatccatgcagctcctatcagcagcaccagggtggcagactcaggactacgatggcaacaagattctgtctcagaga
tgggagctctctctgttgcacagactcagtggtgggtgcatcatagctccctgaggtcttcaccttctgggttcaagtga