

RTKN2-rhotekin 2 Gene

Reference: PTXBC025765

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PLEKHK1; bA531F24.1; rhotekin-2; PH domain-containing family K member 1; pleckstrin homology domain-containing family K member 1		

For research use only.

GENE DESCRIPTION	rhotekin 2
GENE ID	219790
ORIGIN SPECIES	Human
ACCESSIONS	BC025765
VECTOR	pDONR223
NCBI SYMBOL	RTKN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagggccgagcctgaggggtcctgcgtccgctggcggggctcccaccagcaggactgcaacattcaagaaaaatagacttagaattcgaaatcgagaagggaatggaaactcttctctgagcactcagaagatcaagttttacatgcagttaagaatctcatggtgtg
caatgctcgactaatggcctatacatcgagctacagaaattagaagaacagattgcaaatcagactggaagatgtgatgtgaaatttgaagtaagaaacgacagcatgtaaggaaagattgccatcagatattcgaaaccactaatgtggaagactctgatcacttcagcaataaa
gaacgatcacgacgctatgccattttttttattcaaaatgggagcctaattgtgttgatactgatgtgtggaatgtggataaaacaatcacagatatatgttttgaagaatgtaacctattgttaa