

GHRL-ghrelin/obestatin prepropeptide Gene

Reference: PTXBC025791

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MTLRP; appetite-regulating hormone; In2c-preproghrelin; ghrelin, growth hormone secretagogue receptor ligand; ghrelin/obestatin preprohormone; ghrelin/obestatin prepropeptide; growth hormone-releasing peptide; motilin-related peptide; prepro-appetite regu		

For research use only.

GENE DESCRIPTION	ghrelin/obestatin prepropeptide
GENE ID	51738
ORIGIN SPECIES	Human
ACCESSIONS	BC025791
VECTOR	pDONR223
NCBI SYMBOL	GHRL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctccccaggaccgtctgcagctctctgctctcgcatgctctggctggacttggccatggcaggtccagcttctgagccctgaacaccagagagtcagcagagaaaggagtcgaagaagccaccagccaagctgcagccccagctctagcaggtggtccgcccc
gaagatggaggtcaagcagaaggccagaggaatggaagtcgggttcaacgccctttagtgttggaatcaagctgtcaggggttcagtagccagcagcagccagccctgggaagtttcttcaggacatcctctgggaagaggccaaaggccccagccgacaagtga