



Reference: PTXBC025995

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	piwil2-like protein; 80kDa PIWIL2 short isoform; CT80; HILI; PIWIL1L; mili; piwi-like protein 2; Miwi like; cancer/testis antigen 80; piwi-like 2; testicular tissue protein Li 141; piwi like RNA-mediated gene silencing 2		

For research use only.

GENE DESCRIPTION	piwi-like 2 (Drosophila)
GENE ID	55124
ORIGIN SPECIES	Human
ACCESSIONS	BC025995
VECTOR	pDONR223
NCBI SYMBOL	PIWIL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr