

## TMEM217-transmembrane protein 217 Gene

Reference: PTXBC026012

|                               |                                                |                     |                |
|-------------------------------|------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                            | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | C6orf128; dJ355M6.2; transmembrane protein 217 |                     |                |
| <b>For research use only.</b> |                                                |                     |                |

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | transmembrane protein 217                                                                         |
| <b>GENE ID</b>          | 221468                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC026012                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | TMEM217                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaaacagcagcagtggtgtgggatgactgccaaatggcaccgtgtgttcagggtctcaccatcatggccgtagacatgtatctcatctttgaacagaagcacctagggatggcagttgcactgagatcacaccaaagtacaggggtgcaagtaacataataaactcatcatc  
tgcctggagttttaaatcgctctcttctgtcttcatcaccatcctcatcagctgcttctctgtactcagtgatgccagatcttcagggtgctgcatctacattgtctggattttttctatgaaactgcaaacgtcgtatataaatctcaccacaatgacttgacattaaagaggtcagaa  
tcatgcgctggtttggcttgggtgtctctgacagtcagctgttctggatgtcttctgtcatcaactatgccacataacctacaaaaaccggagccagggcaataataatttctacaagagacgaatttctacagcggagattctccacagcagaaataaagattatcaatttcgagtggttca  
gtggctcacacctgggaatccagctacttgagaggcagagggtgctctgtgtaaaacaagtataaaatga