

RAB20-RAB20, member RAS oncogene family Gene

Reference: PTXBC026025

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB20, member RAS oncogene family; ras-related protein Rab-20		
For research use only.			

GENE DESCRIPTION	RAB20, member RAS oncogene family
GENE ID	55647
ORIGIN SPECIES	Human
ACCESSIONS	BC026025
VECTOR	pDONR223
NCBI SYMBOL	RAB20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggaagcccgacagcaagatcgtctctggggacatgaacgtggggaagacgtcgtctgcagcgggtatagagcggcgctccggacacgggtcagcacgggtggcgccgctctacatgaagcagtggtcgtctacaacatctccatctgggacaccgagggcg
ggagcagttccacggcctgggtccatgactgcccggggcgccgccatcctcacctatgatgtgaatcaccggcagagcctggtggagctggaggaccggtctctggcctgacagacacagccagcaagactgcctctcgccatctggggaacaaagtggacctcact
gaggagggggccttggcggccagggagaaggaagagtgagtcaccaatagacgtggggaccgtgtctcccaaggccacctaagcaggtgcagctggaggatgcgtggcccttataaaaaagatcctcaagtacaagatgctggatgagcaggatgtccggccgtgagca
aatgtgctttgagaccagcgccaagaccggctacaatgtggacctctgtttgagacctctttgacctggtggtgccaatgatcttacagcagagagctgagaggccgtcacacacagtgatataatccagtcataagccaccaagaggaccagatctgggtgtgtgcctga