

RBMY2FP-RNA binding motif protein, Y-linked, family 2, member F pseudogene Gene

Reference: PTXBC026077

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	RNA binding motif protein, Y-linked, family 2, member F pseudogene
GENE ID	159162
ORIGIN SPECIES	Human
ACCESSIONS	BC026077
VECTOR	pDONR223
NCBI SYMBOL	RBMY2FP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtagaagcagatcgtctggcaagcatttcgtggcctaataagggaacaatgaaagatgtttaaagcagttatgcgaacatgggtccatcagaagttctttgataaaggatcgaaccagcaatccagaggctttgcatttactctttgagaacactgcagatgctaag
aatgctgccaaacatatgattgaaagcctttggatggaaaagcaataaaagtagaacaagccaagaaaccatctttcaaatgggtgtaggcagagaccaccagcttctcaagaagcagaagccctcaggatgtctgagatcgaaggaagtagtggagggaacaagagatg
gcacccctcacatgaaggacgcttgggtaattgttttaaatgtaagatggaacataggactgaagaaaataagtttgaagatacgaatttctcaattatatttctcttatga