

MOBKL2C-MOB1, Mps One Binder kinase activator-like 2C (yeast) Gene

Reference: PTXBC026078

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MOBKL2C; MOB1E; MOB kinase activator 3C; MOB1, Mps One Binder kinase activator-like 2C; mob1 homolog 2C		

For research use only.

GENE DESCRIPTION	MOB1, Mps One Binder kinase activator-like 2C (yeast)
GENE ID	148932
ORIGIN SPECIES	Human
ACCESSIONS	BC026078
VECTOR	pDONR223
NCBI SYMBOL	MOBKL2C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgtgcctgaagcagggtgtcccaagacaagcgttcggccgggaagcgctttgagccgggcacacagcgctttgagctgtacaagaagcagcagcctctctcaagtcgggcctggacctgcgcagtggtgaggctaccacccggggagaaacatcgacgactgg
atcgccgtgcacgtgggtgactcttcaaccgcatcaacctcatctacggcactatggcgagcgctgcagtgagaccagctgccccgtcatggccggcgggcccccgtacgagtagcctggcaggacgagcgccagtagccggcgccccaagctctctgcgcgcgtatatg
gcgttgcctcatggactggatcgaaggcctcatcaacgacgaagaggtcttccacgcgtgtggagttcccttcccttaagaactccagcaggtctgcaccaagatcctgaccgcctcttccgagctttgtccatgtctacatccaccacttcgatagcatcctcagcatggggcagagg
agcacgtcaacacctgctacaagcacttctacttcatccgcgaggtcagtcgtgtggaccagcgggagctggagccactgaggagatgacagagcggatctgccactga