

C10orf63-chromosome 10 open reading frame 63 Gene

Reference: PTXBC026165

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf63; CFAP106; enkurin; enkurin, TRPC channel interacting protein		
For research use only.			

GENE DESCRIPTION	chromosome 10 open reading frame 63
GENE ID	219670
ORIGIN SPECIES	Human
ACCESSIONS	BC026165
VECTOR	pDONR223
NCBI SYMBOL	C10orf63
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatccaacgtgctctctgagtcatttataacctataccagtgactgaaggagcctcccccagcctctagtgacatatccatttttaaggcaactgtaaaagatgacatgcaaaaagctaaactgcaatgaaactatgggaccagcaaaagttgaagtaccttctccaaaggatttc
ctaaagaacattcaaggaaaaaactctaccaccccaaaaaaactttgatgggaacgtgcccaaaaagcctgctgtgccattgaagactgatcatctgtcatgggaatacagagtggaataaattttataatacaaatgcagctgatatcatatgggagtggtctaaaaagcctaaacca
atttatgttgataaaagaactggagacaagcatgatcttgagccttcaggactagttccaaagtacatcaataaaaggattatgtgtgcacacctgaatacatatgtaagcgaacagaggaaataaaagcccaagaagactatgatctgttatccagaaaaaccttaagaagcagcta
tgaagggtctcctcgatgaagaaaggaggcagttttgcaggggctgaaaaagaactgggaagaggtgcataaagaattccagtcctctcgtctttatagattctataccaagaagatccgcaagcagaggctgggaagaagaatgaacaactagaacacgacattggcataattg
aaaagcacaagattatttatattgccaataacgcatga