

TMEM167A-transmembrane protein 167A Gene

Reference: PTXBC026285

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM167; protein kish-A; transmembrane protein 167; transmembrane protein 167A		
For research use only.			

GENE DESCRIPTION	transmembrane protein 167A
GENE ID	153339
ORIGIN SPECIES	Human
ACCESSIONS	BC026285
VECTOR	pDONR223
NCBI SYMBOL	TMEM167A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgccattttcaatttcagagctattgactgtaactctgtcttatatgtacctgtcttatattgcacccagcctcctggacagaaataaaactggattgttgggtatatttgggaagtgtgccagaattggtgaacggaagagtccttatgttcagtagctgtatagtaatggccttcagcatcctcttcatacagtag