



Reference: PTXBC026309

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD85K; ILT-3; ILT3; LIR-5; LIR5; leukocyte immunoglobulin-like receptor subfamily B member 4; CD85 antigen-like family member K; immunoglobulin-like transcript 3; leucocyte Ig-like receptor B4; leukocyte immunoglobulin-like receptor 5; leukocyte immunoglo		

For research use only.

GENE DESCRIPTION	leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 4
GENE ID	11006
ORIGIN SPECIES	Human
ACCESSIONS	BC026309
VECTOR	pDONR223
NCBI SYMBOL	LILRB4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1