



Reference: PTXBC026314

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Shc3 p51; N-Shc; NSHC; RAI; SHCC; SHC-transforming protein 3; SH2 domain protein C3; SHC (Src homology 2 domain containing) transforming protein 3; SHC-transforming protein C; neuronal Shc; protein Rai; src homology 2 domain-containing transforming protei		

GENE DESCRIPTION	SHC (Src homology 2 domain containing) transforming protein 3
GENE ID	53358
ORIGIN SPECIES	Human
ACCESSIONS	BC026314
VECTOR	pDONR223
NCBI SYMBOL	SHC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]