

ERMN-ermin, ERM-like protein Gene

Reference: PTXBC026345

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA1189; ermin; ermin, ERM-like protein; juxtanodin		

For research use only.

GENE DESCRIPTION	ermin, ERM-like protein
GENE ID	57471
ORIGIN SPECIES	Human
ACCESSIONS	BC026345
VECTOR	pDONR223
NCBI SYMBOL	ERMN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagatgttcgctacattaccaggctgagtgtatggggataaacacacgtgaaacggtcaacaacacactaaatcagtgaggaaatgactgatgtggacagcccccctgccacactacagggtagaacccagctcgaagggtgactcacaaaggaagtcaggaggaa
agaagaaaattacaagggaacatgctgctcaactcatccatggaggacaaaatgctaaaagaaaaccagaagagaaactttattgttcataaggctatcacagatcttctccaagaactagtgctgatgaatgacattcagagaaggccatcagtgaggagaagattcctctgagt
gcagtaaccaggaaataagaagacagaaggagaggattactgagcagcctcacaagagggaagagatgaggacagggaagacaagggtcaccaggcagctgaaattgaatggctgggatttcgaaaacctagccaagctgacatgttaattcctaactgatgaggagcagaagg
tttgggatgaagaaattgatgatgatgatgataattgcaataatgatgaagatgaagttcagtgatagaaatttaagaaaaacatgaagaggtttcctaatttaagagggaagggtgatcaagtgaggactccccactgagcagtgccagttccaagctgtgacacatgatgagcagcc
aaccttagggagaagagtgatatctccagaatgcttattccagatacaatacatcctatcggaatacagaagggaataccaagcaagaattgatgaattcagtgatgatgattataa