



Reference: PTXBC026967

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	palmitoyltransferase ZDHHC5; DHHC5; ZNF375; DHHC-5; membrane-associated DHHC5 zinc finger protein; zinc finger DHHC domain-containing protein 5; zinc finger protein 375; zinc finger, DHHC domain containing 5; zinc finger DHHC-type containing 5		

GENE DESCRIPTION	zinc finger, DHHC-type containing 5
GENE ID	25921
ORIGIN SPECIES	Human
ACCESSIONS	BC026967
VECTOR	pDONR223
NCBI SYMBOL	ZDHHC5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]