

MGC34774-ribosomal protein L13A pseudogene Gene

Reference: PTXBC027852

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	ribosomal protein L13A pseudogene
GENE ID	399670
ORIGIN SPECIES	Human
ACCESSIONS	BC027852
VECTOR	pDONR223
NCBI SYMBOL	MGC34774
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaactcagtcctctgtatcatgatgtcacatccaagccctgcacatcctaggaagattcatcactaagccattactcttcaactaattccatggaaaactcgaactaactcctagcagacatggccgtgtgcacatcaatcatgtcatgaagcaccaggacatcaaggtggctgtgaa
gacagcagaagacagaacctgggtcagtgcccaaggccatcttggggccaccggcagccatcttggccaagtgggctctgctgggaaggaaactgggtgctgagcaccaggggcatcaacttctggaatttctggagaaacaaataaagtatctggccttctctaca
gcagatgagcaccacccctccaaggccgccccgccattcctggggccccagcgcatctttgggtggccttgcagaaccaagtggggccaggctgcctggaccacctatgggtgtttatgggatctcactgcctatga