

THAP10-THAP domain containing 10 Gene

Reference: PTXBC027857

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	THAP domain-containing protein 10; THAP domain containing 10		
For research use only.			

GENE DESCRIPTION	THAP domain containing 10
GENE ID	56906
ORIGIN SPECIES	Human
ACCESSIONS	BC027857
VECTOR	pDONR223
NCBI SYMBOL	THAP10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcccggtgtgtggcgcccgccactgcccgaacaccacacagtcgtgggaagtcgtgtccgtttcccaaggaccggcgctgctgtgggaccgcttcgtgctgggggttgcgcgccgactgttacggaggcaatgaccgctcggtcatctgcttgaccatttgcacca
gcctgtttgacgtctcttcggttatccagaagaacctgcgttccacgcgcctgaggtgtgtggcaggcgccgtgcccacctgcaccgggtgcccccccgccacctaagagggagaggaggagaccagcagccgcttgacacgcgaggagagctccaggcagcca
ggcattctgaggtgccccaggtccagctctgtacacgccccgagctgggaagcaggtgcagcttcacagattacgtgtgaaatgaactgtgcaaacccaacccatgctgataatccatctaatactgtcacttcagttacttacttactgtgaagaaggccagtgataaaagta
cacaaattcttgaagaaggcccgctcaccgtagtgtgggtattcaagccaaagtgaagcggttggaaaaagactgtgtaatgcaactactcagacagagggaattgtgtctagaacttcctctcttgcacatttactccagtgttcagaaacagatacagactgggatatcaagagtgaac
agagtgtttgtcttatgtgctgtacaggtgaaagaagaacatgttaa