

RSPO2-R-spondin 2 homolog (*Xenopus laevis*) Gene

Reference: PTXBC027938

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRISTIN2; R-spondin-2; R-spondin 2 homolog; roof plate-specific spondin-2; R-spondin 2		

For research use only.

GENE DESCRIPTION	R-spondin 2 homolog (<i>Xenopus laevis</i>)
GENE ID	340419
ORIGIN SPECIES	Human
ACCESSIONS	BC027938
VECTOR	pDONR223
NCBI SYMBOL	RSPO2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgccagtagtgagagtgccctgcatctgccatccgggtactatggacacggagcccgatgatgaacagatgtgcaagatgcagaatagaaactgtgattcttcttagcaaaagactttgtaccaagtgcaggtaggctttatttgcatagaggccgttcttgatgaatgtcca
gatggttttgcaccattagaagaacccatggaatgtgtggaaggatgtgaagttgttcattggagcgaatggggaactgttagcagaataatcgacatgttgatttaaatgggtctggaaccagaacacggcaaatgttaaaagccagtgaagacacaataccgtgtccaaccatt
gctgaatccaggagatgcaagatgacaatgaggcattgtccaggagggaagagaacaccaaggcggaaggagaagagggaacaagaaaaagaaagggaagctgatagaaggccaggagcaacacagcgtcttcttagctacagacagagctaaccaataa