

CCL22-chemokine (C-C motif) ligand 22 Gene

Reference: PTXBC027952

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A-152E5.1; ABCD-1; DC/B-CK; MDC; SCYA22; STCP-1; C-C motif chemokine 22; CC chemokine STCP-1; MDC(1-69); chemokine (C-C motif) ligand 22; macrophage-derived chemokine; small inducible cytokine A22; small inducible cytokine subfamily A (Cys-Cys), member 22		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 22
GENE ID	6367
ORIGIN SPECIES	Human
ACCESSIONS	BC027952
VECTOR	pDONR223
NCBI SYMBOL	CCL22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcgctacagactgcactctgtgtgtctctgtctctctgtgtgtggcgcttcaagcaactgaggcagggccctacggcgccaacatggaagacagcgtctgtgtccgtgattacgtccgttaccgtctgccctgcgcgtgtgtgaaacacttctactggacctcagactctgtcccgaggcctggcgtgtgtgtgtaaccttcagggataaggagatctgtgccgatccagagtgccctgggtgaagatgattctcaataagctgagccaatga