



Reference: PTXBC027955

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BDPLT1; BDPLT3; BSS; CD42B; CD42b-alpha; DBPLT3; GPIB; GPIbA; GPIbalpha; VWDP; platelet glycoprotein Ib alpha chain; GP-Ib alpha; antigen CD42b-alpha; glycoprotein Ib (platelet), alpha polypeptide; mutant platelet membrane glycoprotein Ib-alpha; platelet		

GENE DESCRIPTION	glycoprotein Ib (platelet), alpha polypeptide
GENE ID	2811
ORIGIN SPECIES	Human
ACCESSIONS	BC027955
VECTOR	pDONR223
NCBI SYMBOL	GP1BA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1