

GZMH-granzyme H (cathepsin G-like 2, protein h-CCPX) Gene

Reference: PTXBC027974

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCP-X; CGL-2; CSP-C; CTLA1; CTSGL2; cathepsin G-like 2, protein h-CCPX; cytotoxic T-lymphocyte proteinase; cytotoxic T-lymphocyte-associated serine esterase 1; cytotoxin serine protease-C		

For research use only.

GENE DESCRIPTION	granzyme H (cathepsin G-like 2, protein h-CCPX)
GENE ID	2999
ORIGIN SPECIES	Human
ACCESSIONS	BC027974
VECTOR	pDONR223
NCBI SYMBOL	GZMH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagccattctctctgttggcctttcttgacccctgggctgggacagaggatcatcggggccatgaggccaagcccaactcccgcctacatggcctttgttcagtttctgcaagagaagagtcggaagaggtgtggcggcatcctagtggagaaaggactttgtgtgac
agctgtcactgccagggaagctccataaatgtcaccttgggggcccacaatacaaggaacaggagcggacccaagcagtttatccctgtgaaaagaccatccccaatccagcctataatcctaagaacttctccaacgacatcatgctactgagctggagagaagccaagtggac
cacagctgtcggcctctcaggctacatagcagaagccaggtgaagccagggcagctgtgagctgtggctggcgtgggttatgtctcaatgagcacttttagcaaccacactgcagggaagtgtgtgctgacagtgagaaggaactgccagtggtgaactctctccatggcaattacagc
agagccactgagatttgtgtggggatccaaagaacacagaccggttcaagggggactccggggggcccctcgtgtgttaaggacgtagcccaaggtattctctctatggaacaaaaaggacacctccaggagttacatcaaggtctcacacttctgccctggataaagaga
acaatgaagcgctctaa