

C3orf45-chromosome 3 open reading frame 45 Gene

Reference: PTXBC028000

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf45; leucine-rich single-pass membrane protein 2; leucine rich single-pass membrane protein 2		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 45
GENE ID	132228
ORIGIN SPECIES	Human
ACCESSIONS	BC028000
VECTOR	pDONR223
NCBI SYMBOL	C3orf45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccatcattggccccgactgccactgcttgcctgcctgaggagaccaagaagactccgtggcgccaatgatgccagccagaggagcagggggccattggcccccaaccacgtgcatgaggtatgcctgcaccaggtggagtcacagcgaactacatgtggagcagc
acactgcgccctatcactgaagaggcacgacctgggatgagctgctggcgcttttgcgcccgtcactgtgtgccagggctgctgcagccctgtgtacagacgaggagggttcctgctgctcgcgctgctggtgctcactgctagtcgctcactcctgctgtctacctgag
cgtgctgcagagtgaatccctgcgcctcctggcacacacgctccgcacgaggaggagacactactcaaaactccgttggccagcctcagccagcttcggaggctcaactccagtggagcccaagcaccagctga