

RAB39-RAB39, member RAS oncogene family Gene

Reference: PTXBC028064

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB39, member RAS oncogene family; RAB39; ras-related protein Rab-39A; rab-39; rab-related GTP-binding protein; RAB39A, member RAS oncogene family		

For research use only.

GENE DESCRIPTION	RAB39, member RAS oncogene family
GENE ID	54734
ORIGIN SPECIES	Human
ACCESSIONS	BC028064
VECTOR	pDONR223
NCBI SYMBOL	RAB39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaccatctggatctaccagttccgctcatcgtgatcggggactccacgtggcgaagtctgctcctgcaccgcttcaccaggccgcttccccggctgcgctccccgctgcgacccaccgtcggtggacttcttccccgctgctggagatcgagccgggcaag
aggatcaagctacagctctgggacacggcgggacaggagcgggttcagatcaataaccgcatcttattaccgcaactcagttgggtgatttttagtattgacattactaacgcagatctttgaacatgtgaagatfggctagaagaagcaaaatgtatgtacagccatttcggattgtattct
gctagtgggacataaatgtgattagcttcacacgtcaagtacaagggaagaagctgaaaaactgtcagcagactgtggaatgaagtatatagaacctcagcaaaaggatgctacaaatgttgaagaatccttcacaatcttgacgagagacatatgaactattaaaaagggaattt
gtattcaggatggcgggaagggttaaaagtgtttgttccaatactgtgattcttctgaggaagcagtaaaagccaggaaagaatgcttctgctga