

CX3CR1-chemokine (C-X3-C motif) receptor 1 Gene

Reference: PTXBC028078

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCRL1; CMKBRL1; CMKDR1; GPR13; GPRV28; V28; CX3C chemokine receptor 1; C-X3-C CKR-1; CMK-BRL-1; CMK-BRL1; G-protein coupled receptor 13; beta chemokine receptor-like 1; chemokine (C-C) receptor-like 1; chemokine (C-X3-C motif) receptor 1; chemokine (C-X3-		

For research use only.

GENE DESCRIPTION	chemokine (C-X3-C motif) receptor 1
GENE ID	1524
ORIGIN SPECIES	Human
ACCESSIONS	BC028078
VECTOR	pDONR223
NCBI SYMBOL	CX3CR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatcagttccctgaatcagtgacagaaaacttgagtagatgattggctgagggcctgtatattggggacatcgtggtctttgggactgtgttcctgtccatatttactccgtcatcttggccattggcctggtgggaaattgttgtagtgtttgccctaccaacagcaagaagcccaag
agtgacacgacatttacctcctgaacctggccttctgctgctgtttgtgacacattggcctctggaactcactatttgataaatgaaaaggcctccacaatgccatgtgcaaatcactaccgcttctctcatcggctttttggaagcatattctcatcacgctcatcagcattgataggt
acctggccatcgtctgcccgaactcctgaacaacggaccgtgagcatggcgtcaccatcagcctagggcgtctgggcagcagccatttgggtggcagcaccacagttcatgttcacaaagcagaagaaatgaatgccttgggtgactaccccagggtcctcaggaatctggc
ccgtgctcgaatgtggaacaaattttctgctcctactccccctgctcattatgattgtgacttcagaatcagacgctgttttctgcaagaaccacagaagcgaagccattaaactgatcctctggtggtcactcgtgttttctctctgacaccctacaacgttatgatt
tctggagacgcttaagctctatgacttcttcccagttgtgacatgaggaaagctgaggctggcctcagtggtgactgagacgggttcatttagccattgttgcctgaatcctctcatctgatttgcctggggaaggtcagaagataccttaccacgtgatgggaatgctgctgctg
ctgtgtgggcgtcagtcacgttgatttctctcatctgaatcacaaggagcaggcatggaagtgtctgagcagcaatttactaccacacagtgatggagatgcattgctctctctga