



Reference: PTXBC028081

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein C2orf18; C2orf18; ANT2BP; TANGO9; solute carrier family 35 member F6; ANT2-binding protein; transport and golgi organization 9 homolog		

For research use only.

GENE DESCRIPTION	chromosome 2 open reading frame 18
GENE ID	54978
ORIGIN SPECIES	Human
ACCESSIONS	BC028081
VECTOR	pDONR223
NCBI SYMBOL	C2orf18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]